


HP-500™ Plus Vessels

Process up to 14 samples per run with temperatures up to 260 °C
or pressures up to 500 psi.

HP-500 Plus Control Vessel

Need to complete the set	Included with unit		Part number
	0	A. Control Vessel Assembly	431797
1	0	B. Control Spacer	325210
1	0	C. Control Load Disk	430501
1	0	D. Pressure/Temperature Control Cover	431255
1	0	E. Thermowell Nut	325028
1	0	F. Nut, Plug Ferrule	431311
1	0	G. Sapphire Thermowell	326280
1	0	H. Vent Nut	431313
1	0	I. PFA Liner	431110
1	0	J. Vessel Sleeve	431740
0	2	K. Control Screw	430508
0	2	L. Control Frame Module	431755



HP-500 Standard Vessel

Needed to complete the set	Included with set		Part number
		A. Standard Vessel Assembly	431690
5	7	B. Standard Load Disk	430500
12	1	C. Vent Nut	431313
9	4	D. Standard Cover	431250
9	4	E. PFA Liner	431110
1	12	F. Vessel Sleeve	431740
1	12	G. Screw	430507
1	12	H. Standard Frame Module	431755
0	1	Retaining Ring for HP-500 Plus Vessels	431385
0	1	Turntable for HP-500 Plus	431215
0	0	HP-500 Plus Safety Membranes	431300

